Platelet Aggregation Profile

Order Name: PLT AGG
Test Number: 3798203
Revision Date: 10/25/2023

TEST NAME	METHODOLOGY	LOINC CODE
Platelet Aggregation, ADP	Aggregation	5992-3
Platelet Aggregation, Epinephrine	Aggregation	5996-4
Platelet Aggregation, Ristocetin High	Aggregation	5998-0
Platelet Aggregation, Arachadonic Acid	Aggregation	
Platelet Aggregation, Collagen	Aggregation	
Platelet Count for Agglutination	Flow cytometry	777-3
Pathology Report		

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	See Instructions	Whole Blood	Sodium Citrate 3.2% (Blue Top) and EDTA (Lavender Top)	Room Temperature	
Instructions	READ BEFORE COLLECTING SPECIMEN!				
	COLLECTIONS and TESTING MUST BE SCHEDULED with the lab at Ascension St. John Coagulation Department! Please call the coagulation department to make testing arrangements: (918) 744-2500 x15513.				
	If OFF-SITE collection is Authorized by the Coagulations Department, the specimens must reach Ascension St. John Lab within 1 hour of collection and before 1pm. If you cannot arrange for specimens to arrive in this time frame, Do Not Collect Specimen! Patient Must be fasting for at least 8 hours before collection.				
	Collect Both: - Four (2.7mL) 3.2% Sodium (- One (4.5mL) EDTA - Lavende	•			
	Keep specimens whole blood	d (Do Not Spin) - Keep specime	en at Room Temperature!		
	Patient should refrain from asp Greater than 75,000 for accura	71 7	ines or antihistamines for 10 days prior to the test.	Patient should have PLT count	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024

Labcorp Oklahoma, Inc. Test Directory

GENERAL INFORMATION		
Testing Schedule	Mon - Fri (7am - 1pm)	
Expected TAT	1 Day	
Clinical Use	Platelet aggregation studies are done to evaluate platelet function. This is a specialized test and would normally be performed in patients with some indicator of a qualitative platelet disorder.	
CPT Code(s)	85576x5, 8557626	
Lab Section	Coagulation	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2024