

Mycoplasma pneumoniae Antibody

Order Name: **MYCOPLA AB**

Test Number: 5564800

Revision Date: 10/23/2017

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------|-------------------------------|------------|
| Mycoplasma IgG | Indirect Fluorescent Antibody | 6481-6 |
| Mycoplasma IgM | Indirect Fluorescent Antibody | 58708-9 |

SPECIMEN REQUIREMENTS

| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
|---------------------|--|---------------|---------------------------|-----------------------|
| Preferred | 1 mL | Serum | Clot Activator SST | Refrigerated |
| Instructions | Allow specimen to clot completely at room temperature. Separate serum or plasma from cells ASAP or within 2 hours of collection. Stability After separation from cells: Ambient 4hours, Refrigerated 7 days, Frozen 1 months (avoid repeated freeze/thaw cycles). | | | |

GENERAL INFORMATION

| | |
|-------------------------|--|
| Testing Schedule | Mon - Sat |
| Expected TAT | 3 Days |
| Clinical Use | Serologically demonstrates a recent or current infection with M. pneumoniae. Accounts for approximately 20% of all cases of pneumonia. |
| CPT Code(s) | 86738X2 |
| Lab Section | Immunology - Serology |