## Labcorp Oklahoma, Inc. Test Directory

## **Clonazepam Level**

Order Name: CLONAZEPA
Test Number: 4001540
Revision Date: 12/12/2022

TEST NAME		METH	ODOLOGY	LOINC CODE
Clonazepam Level		High F	Performance Liquid Chromatography	3494-2
SPECIMEN REQUIREM	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	4 mL (1.1)	Plasma	EDTA (Lavender Top)	Refrigerated
Alternate 1	4 mL (1.1)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 2	1 mL (0.1)	Plasma	Sodium Heparin (Green Top / No-Gel)	Refrigerated
Instructions	Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gelbarrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.  Specimen Storage: Refrigerate  Specimen Collection: Transfer separated serum or plasma to a plastic transport tube. Collect specimen immediately prior to next dose unless specified otherwise.			

GENERAL INFORMATION		
Expected TAT	2 - 4 days working day.	
Performing Labcorp Test Code	071712	
Notes	Labcorp Test Code: 071712	
CPT Code(s)	80346 (Alt code: G0480)	
Lab Section	Reference Lab	

Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025