Labcorp Oklahoma, Inc. Test Directory

Multiple Sclerosis Panel

icteric CSF specimen

Order Name: MULT SC P
Test Number: 5551560
Revision Date: 10/22/2024

TEST NAME			METHODOLOGY	LOINC CODE
Multiple Sclerosis Panel		Enzyme immunoassay (EIA)		
SPECIMEN REQU	JIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL Serum AND 2.5 mL CSF (See Instructions)	CSF and Serum	Sterile Screwtop Container and Clot Activator SST	Room Temperature
Instructions	Notes: Collect BOTH 2mL (0.8 mL) Serum AND 2.5mL (1.4 mL) CSF Specimen Type: Plastic transport tubes Specimen Storage: Room temperature			
	Specimen Collection: Collect serum AND CSF within 8 hours of each other. Avoid contaminating CSF with blood. Centrifuge CSF at 1500xg for 20 minutes and place the supernatant into separate tubes. Ensure residual fibrin and cellular matter have been removed.			
	Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days			
	Causes for Rejection: Plasma	instead of serum specimen; se	rum and CSF not collected within eight hours of	each other; lipemic, hemolyzed, or

GENERAL INFORMATION			
Expected TAT	4-9 Days		
Clinical Use	Patients with Multiple Sclerosis (MS) have multiple, scarred areas of the brain. Symptoms can initially be mild but typically lead to relapsing or progressive incapacitating neuromotor dysfunction.		
Performing Labcorp Test Code	123390		
CPT Code(s)	82042, 82784x2 82040, 83916, 83873		
Lab Section	Reference Lab		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025