## Labcorp Oklahoma, Inc. Test Directory

## **Gold Level**

Order Name: GOLD
Test Number: 3658375
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
Gold Level			Inductively-Coupled Plasma/Optical Emission Spe	ectroscopy 5659-8
SPECIMEN REQUIRI	EMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	1 mL	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	Specimen Type: Red-top tube OR green-top (heparin) tube  Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.  Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube.  Special Instructions: Testing referred to MedTox Laboratories Inc TC 734			

GENERAL INFORMATION		
Expected TAT	4 - 7 days	
Performing Labcorp Test Code	790750	
Notes	Labcorp Test Code: 790750	
CPT Code(s)	80375	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025