

Selenium Level Urine

Order Name: **SELENIUM U**

Test Number: 3658675

Revision Date: 12/09/2022

TEST NAME	METHODOLOGY	LOINC CODE
Selenium Level Urine	Inductively-Coupled Plasma/Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	20 mL (5 mL)	Urine, Random	Acid Washed, Trace Element Free Container	Room Temperature
Instructions	<p>Specimen Type: Plastic urine container, no preservative</p> <p>Specimen Storage: Maintain specimen at room temperature.</p> <p>Specimen Collection: Optional protocol: Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning). Avoid contact with metal during collection. Screw the lid on securely.</p> <p>Special Instructions: State 24-hour volume on the request form, if applicable. No preservatives. Preservatives used for routine analysis may contain mercuric oxide (ie, Stabilur), which interferes with all metal testing. If both urinalysis and metal testing are ordered, please</p> <p>Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 4 days
Performing Labcorp Test Code	071613
Notes	Labcorp Test Code: 071613
CPT Code(s)	84255
Lab Section	Reference Lab