Transcortin (Corticosteroid-Binding Globulin)

Order Name: **Transcortin**Test Number: 6904943
Revision Date: 12/12/2022

| TEST NAME | | | METHODOLOGY | LOINC CODE |
|---|---|---------------|-------------------------------|-----------------------|
| Transcortin (Corticosteroid-Binding Globulin) | | | Quantitative Radioimmunoassay | |
| SPECIMEN REQUI | REMENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5 mL) | Serum | Clot Activator SST | Frozen |
| Instructions | Notes: 0.5 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: FREEZE Specimen Collection: Separate serum from cells within one hour of collection. Transfer specimen to a plastic transport tube and Freeze. To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED. Special Instructions: Testing is performed at Esoterix Endocrinology Laboratory: ES#500076 | | | |

| GENERAL INFORMATION | | |
|---------------------------------|---------------------------|--|
| Expected TAT | 6 - 10 days | |
| Performing Labcorp Test Code | 500130 | |
| Notes | Labcorp Test Code: 500130 | |
| CPT Code(s) | 84449 | |
| Lab Section | Reference Lab | |
| | | |

Specimen Stability: Ambient: 2 days, Refrigerated : 2 days, Frozen: 300 days

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025