

Immune Complex, C3 Fragments

Order Name: **C3 Fragments**

Test Number: 6907547

Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------------------|--------------------|------------|
| Immune Complex, C3 Fragments | Enzyme Immunoassay | 10864-7 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|---------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.1 mL) | Plasma | EDTA (Lavender Top) | Frozen |
| Alternate 1 | 1 mL (0.1 mL) | Serum | Clot Activator SST | Frozen |
| Instructions | <p>Specimen Type: Red-top tube or gel-barrier tube; lavender-top (EDTA) tube is also acceptable.</p> <p>Specimen Storage: FREEZE at -70(degrees)C or store on dry ice within six hours of collection. Specimen is stable for 14 days if stored at -70 (degrees)C. Storage of specimens at -20(degrees)C will yield erroneous results.</p> <p>Specimen Collection: Transfer specimen to a plastic transport tube before freezing. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p> | | | |

| GENERAL INFORMATION | |
|---------------------|---------------------------|
| Expected TAT | 4 - 8 days |
| Clinical Use | Labcorp Test Code: 123014 |
| CPT Code(s) | 86332 |
| Lab Section | Reference Lab |