## **Cold Agglutinin Antibody**

Order Name: COLD AGG
Test Number: 5500450
Revision Date: 07/10/2025

TEST NAME			METHODOLOGY	LOINC CODE
Cold Agglutinin Antibody			Hemagglutination	5098-9
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5mL)	Serum	Clot Activator (SST or Red N	No-Gel) Refrigerated
Instructions	Specimen: 1mL (0.5mL) Serum from Red-top tube or gel-barrier tube  Collection: Transport blood immediately to laboratory. Incubate at 37°C and allow to clot at 37°C. Separate serum from cells. If there is not a 37°C incubator, allow to clot at room temperature and separate serum from cells immediately after clotting.  Storage Requirements: Refrigerated  Storage Instructions: Do not refrigerate prior to separation of serum from red cells. Store separated serum in refrigerator.  Cause for Rejection: Refrigeration of the specimen before separating serum from cells			

GENERAL INFORMATION			
Expected TAT	3-4 Days		
Clinical Use	Cold agglutinin disease must be considered for all patients with acquired hemolytic anemia who have a positive direct antiglobulin test. Also suspected for patients with chronic liver disease, certain viral infections and Mycoplasma pneumoniae infection.		
Performing Labcorp Test Code	006353		
CPT Code(s)	86157		
Lab Section	Reference Lab		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025