Monospot Test (Mono Test)

Order Name: MONO TEST
Test Number: 5504950
Revision Date: 09/10/2021

| TEST NAME | | METHO | DDOLOGY | LOINC CODE | |
|---------------------------|--|---------------|---------------------------------------|-----------------------|--|
| Monospot Test (Mono Test) | | lmmur | oassay (IA) | 31418-7 | |
| SPECIMEN REQUIREMENTS | | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 1 mL | Serum | Clot Activator SST | Refrigerated | |
| Alternate 1 | 1 mL | Plasma | EDTA (Lavender Top) | Refrigerated | |
| Alternate 2 | 1 mL | Plasma | Lithium Heparin PST (Light Green Top) | Refrigerated | |
| Instructions | Serum should be allowed to clot at room temperature, then centrifuged at 1500 x g for 10 minutes. Serum should be separated ASAP. Refrigerate all specimens at 2°-8°C until ready for testing. If serum or plasma specimens will not be tested within 48hrs, they should be stored at or below -20°C if not tested within 48hrs. Reject criteria: Grossly hemolyzed specimens. | | | | |

| GENERAL INFORMATION | | |
|---------------------|---|--|
| Testing Schedule | Daily | |
| Expected TAT | 1 Day | |
| Clinical Use | Diagnosis of Mononucleosis. Suggest performing an EBV panel if the monospot test is negative. | |
| CPT Code(s) | 86308 | |
| Lab Section | Hematology | |

Note: Further diagnostic workup and/or serology testing is available upon request. Serum is the only acceptable specimen for these serology tests.