

## Calcium Urine Random

Order Name: **CALC R U**  
Test Number: 3002250  
Revision Date: 04/30/2021

TEST NAME	METHODOLOGY	LOINC CODE
Calcium Urine Random	Arsenazo	17862-4

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine specimen. No preservative. Keep refrigerated. 2 - 25 ml of 6 N HCL or 3 - 10 ml of boric acid are acceptable preservatives if collecting with another test. Stability: Ambient 2 days, Refrigerated 4 days, Frozen 3 weeks.			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for identification of abnormal physiologic states causing excess or suppressed excretion of calcium, such as hyperparathyroidism, vitamin D abnormality, diseases that destroy bone, prostate cancer and drug treatment such as thiazide therapy.
CPT Code(s)	82310
Lab Section	Chemistry