

## 11-Deoxycortisol (Compound S)

Order Name: **DEOXYCO11**  
Test Number: 3601100  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
11-Deoxycortisol (Compound S)	Liquid Chromatography/Tandem Mass Spectrometry	1657-6

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated
Alternate 1	1 mL (0.3)	Plasma	EDTA (Lavender Top)	Refrigerated
Alternate 2	1 mL (0.3)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated

**Instructions**

**Specimen Type:** Gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube  
**Specimen Storage:** **Freeze**  
**Specimen Collection:** **Serum/plasma must be separated from cells within 45 minutes of venipuncture, Separate serum/plasma in a plastic transfer tube and Freeze**  
**Special Instructions:** Do not use for Metyrapone test. Testing is performed at Esoterix Endocrinology Laboratory: ES#500132  
**Specimen Stability:** Ambient: 1 day, Refrigerated : 1 day, Frozen: 2 years

### GENERAL INFORMATION

Expected TAT	6 - 15 days
Performing Labcorp Test Code	500171
Notes	Labcorp Test Code: 500171
CPT Code(s)	82634
Lab Section	Reference Lab