Oxalate, Total 24hr Urine

Order Name: OXALATE U
Test Number: 3808300
Revision Date: 12/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Oxalate, Urine - per volume	Quantitative Spectrophotometry	27222-9
Oxalate, Urine - per 24hr	Calculation	2701-1
Creatinine, Urine - per volume	Quantitative Spectrophotometry	20624-3
Creatinine, Urine - per 24hr	Calculation	2162-6

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	4mL (1.5)	Urine, 24-hour	24 hour Urine Container	Frozen	
Instructions	Specimen Type: Plastic LabCorp 24-hour urine container with 30 mL 6N HCl preservative (LabCorp No. 21584). Specimen Storage: Maintain ACIDIFIED urine at room temperature. Acid prevents crystallization and conversion of ascorbate (vitamin C) to oxalate. Specimen Collection: Instruct patient to void at 8 AM (or 8 PM) and discard the specimen. Then collect all urine including the final specimen voide at the end of the 24-hour collection period (ie, 8 AM (or 8 PM) the following day). Label container with patient's name and date and time collection started and finished. Measure and record total urine volume. Mix well. pH must be less than or equal to 3. Preservative not necessary during collection but should be added to urine collection within 24 hours. Special Instructions: State patient's age, sex, and total 24-hour urine volume on the request form and specimen. Specimens submitted without preservative should be acidified after receipt to pH less than or equal to 3 with 6N HCl prior to aliquotting. Specimen Stability: Ambient: 7 days, Refrigerated: 7 days, Frozen: 14 days			onversion of ascorbate (vitamin C) to Il urine including the final specimen voided ent's name and date and time collection Preservative not necessary during	

GENERAL INFORMATION		
Expected TAT	2-5 Days	
Performing Labcorp Test Code	123281	
Notes	Labcorp Test Code: 123281	
CPT Code(s)	83945; 82570	
Lab Section	Reference Lab	