

Chronic Urticaria PD-BATC

Order Name: **CU INDEX**  
Test Number: **5587555**  
Revision Date: **06/03/2025**

TEST NAME	METHODOLOGY	LOINC CODE
Chronic Urticaria PD-BATC	Flow cytometry	75023-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Frozen
Instructions	<b>Specimen:</b> Serum 1 mL (0.5 mL Minimum Volume: 0.5 mL (Note: This volume does not allow for repeat testing.) <b>Container:</b> Gel-barrier tube or red-top tube <b>Collection:</b> <b>Allow to clot then Transfer separated serum to a plastic transport tube and Freeze.</b> <b>Storage Requirements:</b> Frozen <b>Stability Requirements:</b> Room temperature 1 day, Refrigerated 14 days, Frozen 2 months (Freeze/thaw cycles Stable x3) <b>Cause for Rejection:</b> Specimen that are inadequately labeled, lipemic, hemolyzed, cloudy or microbiologically contaminated.			

GENERAL INFORMATION	
Expected TAT	7 - 10 days
Clinical Use	Chronic Urticaria - Pooled Donor Basophil Activation Testing (PD-BAT)
Performing Labcorp Test Code	505800
CPT Code(s)	86849
Lab Section	Reference Lab