

## Deoxycorticosterone

Order Name: **DEOXYSTER**  
Test Number: 3603055  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Deoxycorticosterone	Mass Spectrometry (MS)	1656-8

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Frozen

**Instructions**

**Notes:** 0.5 mL (Note: This volume Does NOT allow for repeat testing.)  
**Specimen Type:** Red-top tube, gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube  
**Specimen Storage:** **FREEZE**  
**Specimen Collection:** **Separate serum from cells within one hour of collection. Transfer specimen to a plastic transport tube and Freeze.** To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.  
**Special Instructions:** Testing is performed at Esoterix Endocrinology Laboratory: ES#500124  
**Specimen Stability:** Ambient: 3 days, Refrigerated : 3 days, Frozen: 2 years

### GENERAL INFORMATION

Expected TAT	6 - 15 days
Performing Labcorp Test Code	500138
Notes	Labcorp Test Code: 500138
CPT Code(s)	82633
Lab Section	Reference Lab