

## Metoprolol Quantitative

Order Name: **Metoprolol Qnt**

Test Number: 3488438

Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Metoprolol Quantitative	Liquid Chromatography/Tandem Mass Spectrometry	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0 mL ( 0.5 mL)	Serum or Plasma	Clot Activator SST	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#126</p>			

### GENERAL INFORMATION

Expected TAT	10-14 days
Clinical Use	Also known as Lopressor
Performing Labcorp Test Code	811547
Notes	Labcorp Test Code: 811547
CPT Code(s)	80299
Lab Section	Reference Lab