

## Mitotane (Lysodren) Level

Order Name: **MITOTANE**

Test Number: 4502125

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Mitotane (Lysodren) Level	Gas chromatography/flame ionization detection (GC/FID)	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated

<b>Instructions</b>	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. <b>Gel-barrier tubes are NOT acceptable.</b></p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc TC 401</p>			
---------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--	--	--

### GENERAL INFORMATION

Expected TAT	7 - 10 days
Performing Labcorp Test Code	808159
Notes	Labcorp Test Code: 808159
CPT Code(s)	80299
Lab Section	Reference Lab