

Pyruvic Acid, CSF

Order Name: **PYRU CSF**

Test Number: 6905131

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Pyruvic Acid, CSF	Quantitative Enzymatic	14122-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	CSF (Cerebrospinal Fluid)	Sterile Screwtop Container	Frozen
Instructions	<p>Specimen Type: Sterile screw capped vial</p> <p>Specimen Storage: Freeze as soon as possible after collection</p> <p>Specimen Collection: Collect entire sample into a single sterile tube. Immediately after collection, add a 2:1 ratio (acid:sample) of COLD 8% perchloric acid. Vortex for 30 seconds, and place on ice for 5 minutes. Store frozen at -20 degrees Celsius and ship frozen.</p> <p>Special Instructions: Testing referred to Medical Neurogenetic LLC MNEGA#MET11</p> <p>Specimen Stability: Ambient: unstable (stability provided by manufacturer or li, Refrigerated : 24 hours (stability provided by manufacturer or li, Frozen: -20 C = 72 hours; -80 C = Indefinite (stability pr</p>			

GENERAL INFORMATION	
Expected TAT	10 - 14 days
Performing Labcorp Test Code	620045
Notes	Labcorp Test Code: 620045
CPT Code(s)	84210
Lab Section	Reference Lab