

Diabetes Group

Order Name: **Diabetes Group**
Test Number: 2019200
Revision Date: 05/03/2018

| TEST NAME | METHODOLOGY | LOINC CODE |
|--|-------------|------------|
| Glucose | Hexokinase | 2345-7 |
| Hemoglobin A1C (Glycosylated) with Mean Glucose | | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|-------------------------|-------------------------|-------------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | See Instructions | See Instructions | See Instructions | See Instructions |
| Instructions | <p>Glucose Level: Preferred 1 mL (0.5) Plasma Sodium fluoride (gray top) or 1 mL (0.5) Plasma Lithium Heparin PST (Light Green Top) or 1 mL (0.5) Serum Clot Activator SST - Transport/Storage: Refrigerated. Lithium heparin and clot tube must be spun within 2 hours for the integrity of the specimen. Stability: Ambient: 8 hours; Refrigerated: 72 hours Fluoride preserved plasma samples are stable at room temperature for 24 hours.</p> <hr/> <p>Hemoglobin A1c w/ Mean Glucose: 1 mL (0.5 mL) Whole Blood EDTA (Lavender Top) Refrigerated Patient does not need to be fasting. Specimen stability: Ambient 8 hours. Refrigerated 7 days, Frozen Unacceptable.</p> | | | |

| GENERAL INFORMATION | |
|-------------------------|-----------------------|
| Testing Schedule | Sun-Sat |
| Expected TAT | 1-2 Days |
| CPT Code(s) | 82947, 83036 |
| Lab Section | See Individual Assays |