## Labcorp Oklahoma, Inc. Test Directory

## Parvovirus B19 DNA, PCR, Quantitative

Order Name: **PARVO QT** Test Number: 3613445 Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE	
Parvovirus B-19 B	y PCR		Polymerase Chain Reaction	9572-9	
Parvovirus B-19 PCR log10copy/mL		Polymerase Chain Reaction	N/A		
SPECIMEN REQUIREMENTS					
Crassing an		Cracing and Turne	Creating Container	Transmert Freisingen ant	

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Frozen
Alternate 1	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Alternate 2	2 mL (0.5 mL)	Serum	Clot Activator SST	Frozen
Instructions	Specimen Type: Lavender-top (EDTA) tube, red-top tube OR gel-barrier tube			

Specimen Storage: FREEZE(preferred) OR Refrigerated

**Specimen Collection:** Collect specimen in lavender-top (EDTA) tube, red-top tube or gel-barrier tube. Centrifuge specimen within 24 hours of collection. Remove serum or plasma, transfer specimen to polypropylene screw-cap tube, and FREEZE. Ship frozen. To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED. **Specimen Stability:** Ambient: Not Available, Refrigerated : 7 days, Frozen: 90 days

Prompt code: [4182051] - Specimen Source

GENERAL INFORMATION			
Expected TAT	3 - 5 days		
Clinical Use	Detection and quantitation of parvovirus in plasma or serum specimens as an aid in the diagnosis and management of parvovirus infections		
Performing Labcorp Test Code	139326		
Notes	Labcorp Test Code: 139326		
CPT Code(s)	87799		
Lab Section	Reference Lab		