

Urea Nitrogen, Blood Post Dialysis

Order Name: **BUN POST**

Test Number: 2009750

Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Urea Nitrogen, Blood Post Dialysis	Urease/GLDH	11064-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Stability: Ambient 7 days, Refrigerated 7 days, Frozen 1 year			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Post dialysis. Useful in the monitoring of dialysis.
CPT Code(s)	84520
Lab Section	Chemistry