## Prostate Specific Antigen (PSA) Total

Order Name: PSA EQ
Test Number: 2012225
Revision Date: 10/28/2025

TEST NAME			METHODOLOGY	LOINC CODE	
Prostate Specific Antigen (PSA) Total			Chemiluminescence Assay	2857-1	
SPECIMEN REQU	JIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated	
Instructions	Remove Serum from cells/gel and freeze if testing will be delayed >24hrs (or w/in 3 hours if testing PSA Free)				
	Stability: Ambient 1 day. Refri	gerated 7 days. Frozen > 7 days.	ays.		

GENERAL INFORMATION			
Testing Schedule	Mon - Fri		
Expected TAT	1-3 days		
Clinical Use	Useful for screening for prostate cancer. Useful for monitoring patients with a history of prostate cancer.		
Notes	Testing performed via Abbott Architect Immunoassay method. Results obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy.		
CPT Code(s)	84153		
Lab Section	Chemistry		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025