

Cyclosporine Level Peak

Order Name: **CYCLO PEAK**

Test Number: 2015325

Revision Date: 05/03/2021

| TEST NAME               | METHODOLOGY                                | LOINC CODE |
|-------------------------|--|------------|
| Cyclosporine Level Peak | Chemiluminescent Micropartical Immunoassay | 53834-8    |

| SPECIMEN REQUIREMENTS |  |               |                     |                       |
|-----------------------|--|---------------|---------------------|-----------------------|
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container  | Transport Environment |
| Preferred             | 2 mL (1.0)   | Whole Blood   | EDTA (Lavender Top) | Refrigerated          |
| Instructions          | Peak: draw 2 hours after medication.<br>Stability: Ambient 24 hours. Refrigerated 7 days. Frozen > 7 days. |               |                     |                       |

| GENERAL INFORMATION |   |
|---------------------|---|
| Testing Schedule    | Mon - Fri   |
| Expected TAT        | 1-2 days  |
| Clinical Use        | This test is used to monitor blood concentration of Cyclosporine A. |
| CPT Code(s)         | 80158   |
| Lab Section         | Chemistry   |