Labcorp Oklahoma, Inc. Test Directory

CA 125, Serous Fluid

Order Name: SRS CA125
Test Number: 2015430
Revision Date: 10/28/2025

TEST NAME			METHODOLOGY	LOINC CODE
CA 125, Serous Fluid			Chemiluminescence Assay	11210-2
SPECIMEN REQU	JIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serous Fluid	Sterile Screwtop Container	Refrigerated
Instructions	Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Stability: Ambient 1 day. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION		
Testing Schedule	Mon - Fri	
Expected TAT	1-3 days	
Clinical Use	Useful in the follow-up management of patients undergoing cancer therapy, especially for ovarian carcinoma.	
Notes	Testing performed via Abbott Alinity Immunoassay method. Results obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy.	
CPT Code(s)	86304	
Lab Section	Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025