## CA 125, Serous Fluid

Order Name: SRS CA125 Test Number: 2015430 Revision Date: 06/14/2018

TEST NAME			METHODOLOGY	LOINC CODE
CA 125, Serous Fluid			Chemiluminescence Assay	11210-2
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serous Fluid	Sterile Screwtop Container	Refrigerated
Instructions	ons Venous blood is often drawn simultaneously. Note fluid type on requisition and container.			
	Stability: Ambient 1 day. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION		
Testing Schedule	Mon - Fri	
Expected TAT	1-3 days	
Clinical Use	Useful in the follow-up management of patients undergoing cancer therapy, especially for ovarian carcinoma.	
CPT Code(s)	86304	
Lab Section	Chemistry	