CA 19-9 (Carbohydrate Antigen 19-9)

Order Name: CA 19-9
Test Number: 2024050
Revision Date: 10/28/2025

TEST NAME			METHODOLOGY	LOINC CODE
CA 19-9 (Carbohydrate Antigen 19-9)			Chemiluminescence Assay	24108-3
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2mL (1mL)	Serum	Clot Activator SST	Refrigerated
Alternate 1	2mL (1mL)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Stability: Ambient 1 day. Refrigerated 7 days. Frozen > 7 days. Unacceptable: gross hemolysis			

GENERAL INFORMATION			
Testing Schedule	Mon-Fri		
Expected TAT	1-3 days		
Clinical Use	Potentially useful adjunct for diagnosis and monitoring of pancreatic cancer.		
Notes	Testing performed via Abbott Alinity Immunoassay method. Results obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy.		
CPT Code(s)	86301		
Lab Section	Chemistry		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025