

Potassium Urine Random

Order Name: **POTAS R U**
Test Number: 3001350
Revision Date: 05/28/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------------|-------------------------|------------|
| Potassium Urine Random | Ion-Selective Electrode | 2828-2 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|-------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Urine, Random | Sterile Urine container | Refrigerated |
| Instructions | Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 24 hours, Refrigerated 2 months, Frozen 1 year | | | |

| GENERAL INFORMATION | |
|---------------------|---|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful in the management of kidney disease. |
| CPT Code(s) | 84133 |
| Lab Section | Chemistry |