## Protein Body Fluid (Synovial/Serous)

Order Name: SRS PROT
Test Number: 3500450
Revision Date: 06/15/2017

TEST NAME			METHODOLOGY	LOINC CODE
Protein Body Fluid (Synovial/Serous)		Biuret	2881-1	
SPECIMEN REQUIRE	EMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Body Fluid	Sterile Screwtop Container	Refrigerated
Alternate 1	1 mL (0.5)	Body Fluid	Lithium Heparin (Dark Green Top No-Gel)	/ Refrigerated
Instructions	Place body fluid in a Sterile Screwtop container or Lithium Heparin Green top tube with No-Gel.  Venous blood is often drawn simultaneously. <b>Note fluid type on requisition and container</b> .  Specimen stability: Ambient 8 hours. Refrigerated 7 days.			

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1-2 days	
Clinical Use	Used in evaluating protein nutritional status and protein altering diseases.	
CPT Code(s)	84157	
Lab Section	Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025