Triglyceride Serous Fluid

Order Name: SRS TRIG
Test Number: 3502150
Revision Date: 06/12/2017

TEST NAME			METHODOLOGY	LOINC CODE	
Triglyceride Serous Fluid		Glycerol Phosphate Oxidase	12228-3		
SPECIMEN REQUIRI	EMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.5 mL)	Serous Fluid	Sterile Screwtop Container	Refrigerated	
Instructions	Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen stability: Ambient 8 hours. Refrigerated 7 days.				
GENERAL INFORMATION					
Testing Schedule	Daily				

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1-2 days	
Clinical Use	Useful for evaluation of risk for cardiovascular disease and in evaluation of hyperlipidemia.	
CPT Code(s)	84478	
Lab Section	Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025