## Follicle Stimulating Hormone - FSH, Serum

Order Name: FSH

Test Number: 3601200 Revision Date: 05/20/2021

TEST NAME			METHODOLOGY	LOINC CODE
Follicle Stimulating Hormone - FSH, Serum			Chemiluminescence Assay	15067-2
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Remove Serum/Plasma from cells/gel if testing delayed >24hrs.  Stability: Ambient 1 day. Refrigerated 7 days. Frozen 12 months.			

GENERAL INFORMATION		
Testing Schedule	Mon - Fri	
Expected TAT	1-2 days	
Clinical Use	Useful for evaluation of menstrual irregularities, work-up of patients with suspected hypogonadism, prediction of ovulation, evaluation of infertility and the diagnosis of pituitary disorders.	
CPT Code(s)	83001	
Lab Section	Chemistry	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025