Viscosity, Serum

Order Name: VISCOSITY
Test Number: 2005770
Revision Date: 12/12/2022

| TEST NAME | | | METHODOLOGY | LOINC CODE |
|-----------------------|--|---------------|--------------------|-----------------------|
| Viscosity, Serum | Cone-Plate Viscometer | | | |
| SPECIMEN REQUIREMENTS | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Serum | Clot Activator SST | Refrigerated |
| Instructions | Notes: 0.5 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Refrigerate Specimen Collection: Not Available Specimen Stability: Ambient: 14 days, Refrigerated: 14 days, Frozen: 14 days | | | |

| GENERAL INFORMATION | | |
|---------------------------------|---------------------------|--|
| Expected TAT | 3 - 5 days shifts. | |
| Performing Labcorp Test Code | 004861 | |
| Notes | Labcorp Test Code: 004861 | |
| CPT Code(s) | 85810 | |
| Lab Section | Reference Lab | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025