

Chlorpromazine (Thorazine) Level

Order Name: **THORAZINE**

Test Number: 2008650

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Chlorpromazine (Thorazine) Level	Liquid Chromatography/Tandem Mass Spectrometry	3471-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Alternate 2	2 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Instructions	<p>Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.</p> <p>Specimen Storage: Room Temperature</p> <p>Specimen Collection: Transfer separated serum or plasma to a plastic transport tube.</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 5 days
Performing Labcorp Test Code	072132
Notes	Labcorp Test Code: 072132
CPT Code(s)	80342
Lab Section	Reference Lab