Beta-2, Transferrin BF

Preferred

Instructions

Order Name: BETA-2 TRA
Test Number: 3656675
Revision Date: 10/30/2025

Refrigerated

Sterile Screwtop Container

TEST NAME			METHODOLOGY	LOINC CODE
Beta-2, Transferrin BF			Immunofixation Electrophoresis	
SPECIMEN REQU	UIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment

Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens!

Body Fluid

Notes: 1 mL (Note: This volume Does NOT allow for repeat testing).

Specimen Type: Sterile container Specimen Storage: Refrigerated

2 mL

Special Instructions: Testing referred to ARUP TC 50047

Specimen Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks

GENERAL INFORMATION		
Expected TAT	4 - 8 days	
Clinical Use	This test is to detect CSF in body fluids such as sinus or ear drainage.	
Performing Labcorp Test Code	829030	
CPT Code(s)	86335	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025