## Cryofibrinogen

Order Name: CRYOFIBRIN
Test Number: 5221675
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
Cryofibrinogen			Precipitation	
SPECIMEN REQUIR	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (1mL)	Plasma	Sodium Citrate 3.2% (Blue Top)	See Instructions
Instructions	Specimen Type: Blue-top (sodium citrate) tube, lavender-top (EDTA) tube, or gray-top (sodium fluoride) tube  Specimen Storage: Maintain specimen at room temperature.  Specimen Collection: Specimen must be drawn in a prewarmed tube and kept at 37(degrees)C. Separate plasma from cells immediately and transfer the plasma into a clean transport tube. LABEL TRANSPORT TUBE "PREWARMED."  Specimen Stability: Ambient: Not Available. Refrigerated: Not Available. Frozen: Not Available			

GENERAL INFORMATION		
Expected TAT	3 - 5 days later.	
Performing Labcorp Test Code	080317	
Notes	Labcorp Test Code: 080317	
CPT Code(s)	82585	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025