Labcorp Oklahoma, Inc. Test Directory

HazeInut Meat IgG

Preferred

Instructions

Order Name: HAZELNUT G

Room Temperature

Test Number: 5500463

Revision Date: 12/14/2022

TEST NAME			METHODOLOGY	LOINC CODE
Hazelnut Meat IgG			Enzyme immunoassay (EIA)	
SPECIMEN REQ	UIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment

Specimen Type: Red-top tube or gel-barrier tube

Specimen Storage: Room temperature **Specimen Collection:** Not Available

0.3 mL (0.1 mL)

Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only

Clot Activator SST

one assay method be used consistently to monitor each patient's course of therapy.

Serum

Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available

GENERAL INFORMATION			
Expected TAT	2 - 4 days		
Clinical Use	The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test. The reference range listed on the report is the lower limit of quantitation for the assay.		
Performing Labcorp Test 605043 Code			
Notes	Labcorp Test Code: 605043		
CPT Code(s)	86001		
Lab Section	Reference Lab		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025