

Butalbital Level (Fiorinal)

Order Name: Butalbital

Test Number: 3608250

Revision Date: 11/13/2024

TEST NAME	METHODOLOGY	LOINC CODE
Butalbital Level (Fiorinal)	Liquid Chromatography/Tandem Mass Spectrometry	6895-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0mL (0.5mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2.0mL (0.5mL)	Plasma	Sodium Heparin (Green Top, No-Gel)	Refrigerated
Instructions	<p>Specimen: 2mL (0.5mL) Serum or Plasma - Refrigerated</p> <p>Container: Red-Top (no additive) Tube or Green-Top (heparin) Tube. Gel- barrier tubes are not recommended.</p> <p>Testing referred to MedTox Laboratories Inc. MX#22</p> <p>Collection: Transfer separated serum or plasma to a plastic transport tube.</p> <p>Storage Instructions: Submission/transport (<3 days) Room temperature. For storage beyond three days, specimen should be Refrigerated or Frozen.</p>			

GENERAL INFORMATION	
Expected TAT	Set up Monday thru Friday on 1st shift; reported 1-2 Days from setup
Clinical Use	Also known as: Axocet, Lotusate, and Sandoptal
Performing Labcorp Test Code	716320
CPT Code(s)	80345
Lab Section	Reference Lab