Labcorp Oklahoma, Inc. Test Directory

Methaqualone Quantitative, Urine

Order Name: Methaqual Ur
Test Number: 6900779
Revision Date: 04/24/2021

TEST NAME			METHODOLOGY	LOINC CODE	
Methaqualone Quantitative, Urine			Gas Chromatography/Mass Spectrometry (Go	C/MS) 3787-9	
SPECIMEN REQUIRE	EMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	20 mL	Urine, Random	Sterile Screwtop Container	Room Temperature	
Instructions	Specimen Type: Use plastic urine drug bottle and tamper-evident seal for forensic specimen. Collection kits are available by request from the laboratory. Specimen Storage: Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate. Specimen Collection: Urine temperature monitoring is recommended for samples to be tested for medicolegal purposes.				
	Special Instructions: Chain-of-custody documentation is required for samples submitted for pre-employment, random employee testing, and forensic purposes. For other applications, use the standard request form. Please mark chain-of-custody test number on the request form. Specimen Stability: Ambient: Not Available, Refrigerated: Not Available, Frozen: Not Available				

GENERAL INFORMATION	
Performing Labcorp Test Code	071332
Notes	Labcorp Test Code: 071332
CPT Code(s)	80375
Lab Section	Reference Lab

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025