

Deoxypyridinoline Crosslinks (Pyrilinks-D)

Order Name: **Pyrilinks-D**

Test Number: 6904991

Revision Date: 04/10/2017

TEST NAME	METHODOLOGY	LOINC CODE
Deoxypyridinoline Crosslinks (Pyrilinks-D)	Immunochemiluminometric (ICMA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, First Void	Sterile Urine container	Frozen
Instructions	<p>Specimen Type: Plastic urine container</p> <p>Specimen Storage: Refrigerate</p> <p>Specimen Collection: Collection before 10 AM is recommended (eg, first morning void). Do NOT add preservative.</p> <p>Special Instructions: The account must submit the patient's Social Security number to monitor.</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	4-5 days
Performing Labcorp Test Code	511105
Notes	Labcorp Test Code: 511105
CPT Code(s)	82570, 82523
Lab Section	Reference Lab