

Periprosthetic Joint Infection (PJI) Detection (Synovasure or Alpha Defensin)

Order Name: **SYNOVA PJI**
Test Number: 6905105
Revision Date: 03/31/2022

TEST NAME	METHODOLOGY	LOINC CODE
PJI Source	Prompt	31208-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Synovial Fluid	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Collect Synovial fluid in Plain Red No-Gel Clot tube. Transport 1mL(0.5mL)synovial fluid in the original collection tube. Separate specimens must be submitted when multiple tests are ordered. Stability: Ambient: 72 hours; Refrigerated: 4 days; Frozen: Unacceptable Due to limited stability, only collect Mon-Wed and not the day before a holiday			

GENERAL INFORMATION	
Testing Schedule	Varies
Expected TAT	3-5 days
Notes	Reference Lab: ARUP Test Code: 2013008 Click Here to view information on the ARUP website.
CPT Code(s)	83516, 86140
Lab Section	Reference Lab