

Flunitrazepam (Rohypnol) Scr.U

Order Name: **FNZ Scr Ur**

Test Number: 5194624

Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|------------------|------------|
| Flunitrazepam (Rohypnol) Scr.U | Immunoassay (IA) | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|-------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 10 mL (3 mL) | Urine, Random | Sterile Urine container | Room Temperature |
| Instructions | <p>Notes: 3 mL (Note: This volume Does NOT allow for repeat testing).</p> <p>Specimen Type: Plastic urine container without preservative</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Special Instructions: Testing referred to MedTox Laboratories Inc MX#436</p> | | | |

| GENERAL INFORMATION | |
|------------------------------|---------------------------|
| Expected TAT | 3-10 Days |
| Performing Labcorp Test Code | 808417 |
| Notes | Labcorp Test Code: 808417 |
| CPT Code(s) | 80307 |
| Lab Section | Reference Lab |