

Cat Scratch Disease Antibody (Bartonella)

Order Name: **CATSCRATCH**

Test Number: 5590000

Revision Date: 12/26/2017

TEST NAME	METHODOLOGY	LOINC CODE
Bartonella henselae Antibody IgM	Indirect Fluorescent Antibody	32576-1
Bartonella henselae Antibody IgG	Indirect Fluorescent Antibody	32575-3
Bartonella quintana Antibody IgG	Indirect Fluorescent Antibody	9360-9
Bartonella quintana Antibody IgM	Indirect Fluorescent Antibody	32578-7

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (1)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen should be collected in a Clot Activator SST with a gel barrier. Ambient specimens are unacceptable.			

GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	3 Days
Clinical Use	Positive serology to Bartonella henselae and Bartonella quintana assist in the diagnosis of Cat Scratch disease.
CPT Code(s)	86611X4
Lab Section	Immunology - Serology