

Cyclosporine Level Peak

Order Name: **CYCLO PEAK**

Test Number: 2015325

Revision Date: 05/03/2021

TEST NAME	METHODOLOGY	LOINC CODE
Cyclosporine Level Peak	Chemiluminescent Micropartical Immunoassay	53834-8

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1.0)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Instructions	Peak: draw 2 hours after medication. Stability: Ambient 24 hours. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION

Testing Schedule	Mon - Fri
Expected TAT	1-2 days
Clinical Use	This test is used to monitor blood concentration of Cyclosporine A.
CPT Code(s)	80158
Lab Section	Chemistry