

pH Serous Fluid

Order Name: **SRS PH**

Test Number: 3500350

Revision Date: 03/06/2020

TEST NAME	METHODOLOGY	LOINC CODE
pH Serous Fluid		2748-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serous Fluid	Sterile Screwtop Container	On Ice
Instructions	Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen must be on ice after collection. Deliver to lab immediately.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful in assessing acid-base balance.
CPT Code(s)	83986
Lab Section	Chemistry