

## Estriol, Serum

Order Name: **ESTRIOL**  
Test Number: 3801500  
Revision Date: 12/12/2022

| TEST NAME      | METHODOLOGY                    | LOINC CODE |
|----------------|--------------------------------|------------|
| Estriol, Serum | Immunochemiluminometric (ICMA) | 2251-7     |

### SPECIMEN REQUIREMENTS

| Specimen            | Specimen Volume (min)   | Specimen Type | Specimen Container | Transport Environment |
|---------------------|---|---------------|--------------------|-----------------------|
| Preferred           | 1 mL (0.3)  | Serum         | Clot Activator SST | Refrigerated          |
| <b>Instructions</b> | <p><b>Notes:</b> 0.3 mL (Note: This volume Does NOT allow for repeat testing.)<br/> <b>Specimen Type:</b> Red-top tube or gel-barrier tube<br/> <b>Specimen Storage:</b> Refrigerate<br/> <b>Specimen Collection:</b> If red-top tube is used, transfer separated serum to a plastic transport tube. All specimens should be drawn at the same time of day to allow comparison of values.<br/> <b>Special Instructions:</b> State weeks of gestation on the request form. Usually done as serial measurement.<br/> <b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p> |               |                    |                       |

### GENERAL INFORMATION

|                              |                           |
|------------------------------|---------------------------|
| Expected TAT                 | 2 - 4 days                |
| Performing Labcorp Test Code | 004614                    |
| Notes                        | Labcorp Test Code: 004614 |
| CPT Code(s)                  | 82677                     |
| Lab Section                  | Reference Lab             |