

Amitriptyline

Order Name: **AMITRIPTL**

Test Number: 4302455

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Amitriptyline		3333-2
Nortriptyline (Aventyl)	Quant Liquid Chromatography-Tandem Mass Spectrometry	3872-9
Amitriptyline/Nortriptyline		3335-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	1 mL (0.3 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Room Temperature
Alternate 2	1 mL (0.3 mL)	Plasma	EDTA (Lavender Top)	Room Temperature
Instructions	Notes: 0.3 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Specimen Storage: Room Temperature Specimen Collection: Transfer separated serum or plasma to a plastic transport tube. For therapeutic monitoring, collect specimen immediately prior to next dose unless specified otherwise. Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days			

GENERAL INFORMATION	
Expected TAT	2 - 5 days Sunday on 3rd shift; reported within one working day.
Performing Labcorp Test Code	007476
Notes	Labcorp Test Code: 007476
CPT Code(s)	80335
Lab Section	Reference Lab