

Vaginosis Profile Saline

Order Name: **VAG PROF**
Test Number: 2915425
Revision Date: 09/17/2018

TEST NAME	METHODOLOGY	LOINC CODE
Ph of vaginal discharge		48037-6
Whiff test	Amine	46730-8
Gram stain		
Trichomonas Microscopic Examination	Microscopy	32766-8
Vaginal Yeast Examination	Microscopy	32765-0
Clue Cell Examination	Microscopy	32764-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	Swab	Sterile Saline	Room Temperature
Instructions	<p>Obtain specimen from mucosal surface and place swab in saline tube. Samples contaminated with preparations containing iodine or by the immediate prior use of vaginal lubricants are not recommended. For optimal results, testing should be performed as soon as possible following collection; specimens must be <24 hours old upon arrival.</p> <p>Note: The Vaginal Profile from Swab (with Trichomonas Antigen), is recommended for specimens being testing >2 hours after collection. Please refer to test 2915445 for collection information.</p>			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Determines bacterial vaginosis or yeast vaginitis
CPT Code(s)	83986; 82120; 87210; 87205
Lab Section	Microbiology