## Cow Dander Allergen

Order Name: COW DAND Test Number: 5609150 Revision Date: 02/11/2013

| TEST NAME   |  |               | METHODOLOGY        | LOINC CODE              |
|---|--|---------------|--------------------|-------------------------|
| Cow Dander Allergen   |  | ImmunoCAP     | 15661-2            |                         |
| SPECIMEN REQUIREMENTS   |  |               |                    |                         |
| Specimen  | Specimen Volume (min)  | Specimen Type | Specimen Container | Transport Environment   |
| Preferred   | 0.4 mL (0.2 mL)  | Serum         | Clot Activator SST | Room Temperature        |
| Instructions  | tions Specimen Type: Red-top tube or gel-barrier tube,         |               |                    |                         |
|   | Separate Serum from Cells into a screwtop transport container. |               |                    |                         |
| Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stal |  |               |                    | /thaw cycles Stable x3) |

| GENERAL INFORMATION |               |  |
|---------------------|---------------|--|
| Expected TAT        | 3-5 days      |  |
| CPT Code(s)         | 86003         |  |
| Lab Section         | Reference Lab |  |