

Animal and Feather Allergen Panel

Order Name: **A ANIML PN**

Test Number: **5616750**

Revision Date: **08/01/2023**

TEST NAME	METHODOLOGY	LOINC CODE
Cat Dander Allergen	ImmunoCAP	15609-1
Dog Dander Allergen	ImmunoCAP	15685-1
Cow Dander Allergen	ImmunoCAP	15661-2
Horse Dander Allergen	ImmunoCAP	15784-2
Feathers Allergen Panel	ImmunoCAP	31161-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<b>Specimen Type:</b> Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. <b>Stability Requirements:</b> Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 Days
Notes	Performed at Labcorp Phoenix as individual allergens
CPT Code(s)	86003x5
Lab Section	Reference Lab