

Butalbital Level (Fiorinal)

Order Name: **Butalbital**
Test Number: 3608250
Revision Date: 11/13/2024

TEST NAME	METHODOLOGY	LOINC CODE
Butalbital Level (Fiorinal)	Liquid Chromatography/Tandem Mass Spectrometry	6895-7

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0mL (0.5mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2.0mL (0.5mL)	Plasma	Sodium Heparin (Green Top, No-Gel)	Refrigerated
Instructions	Specimen: 2mL (0.5mL) Serum or Plasma - Refrigerated Container: Red-Top (no additive) Tube or Green-Top (heparin) Tube. Gel- barrier tubes are not recommended. Testing referred to MedTox Laboratories Inc. MX#22 Collection: Transfer separated serum or plasma to a plastic transport tube. Storage Instructions: Submission/transport (<3 days) Room temperature. For storage beyond three days, specimen should be Refrigerated or Frozen.			

GENERAL INFORMATION

Expected TAT	Set up Monday thru Friday on 1st shift; reported 1-2 Days from setup
Clinical Use	Also known as: Axocet, Lotusate, and Sandoptal
Performing Labcorp Test Code	716320
CPT Code(s)	80345
Lab Section	Reference Lab