Labcorp Oklahoma, Inc. Test Directory

Toluene, Blood

Order Name: **TOLUENE**Test Number: 3610275
Revision Date: 12/12/2022

LOINC

IESI NAME			METHODOLOGY	CODE
Toluene, Blood			Headspace Gas Chromatography/Flame Ioniza FID)	ation Detection (GC- 9637-0
SPECIMEN REQUI	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	20 mL (7 mL)	Whole Blood	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	Specimen Type: Green-top (heparin) tube Specimen Storage: Refrigerate; specimen may also be frozen. Specimen Collection: Collect adequate number of tubes to obtain the required specimen volume. Tubes should be completely filled to prevent loss of volatile compound into headspace. Collect prior to last shift of workweek.			

GENERAL INFORMATION		
Expected TAT	10 - 14 days	
Performing Labcorp Test Code	790580	
Notes	Labcorp Test Code: 790580	
CPT Code(s)	84600	
Lab Section	Reference Lab	

Special Instructions: Testing referred to MedTox Laboratories Inc TC 70123

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025